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The HTRV-A split tube furnaces have a maximum operating temperature of 1700 °C.

The split heating module allows either easy positioning of the work tube or positioning around reactors which have fixed end flanges. The split design may also allow faster cooling of samples. The control thermocouple is fitted in the centre of the heating zone. Cooling channels are engineered into the housing to aid with convection cooling of the outer case.

The two furnace chamber halves consist of high grade insulation plates with vertically hanging MoSi₂ heating elements. A safety switch protects the operator by switching off the heating elements once the furnace is opened.

The tube furnace is supplied without a stand, allowing customers to build them into their own equipment. Optional 'L' stands are available allowing the furnaces to be self supporting.

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carbonizare, creștere cristale, călire, degazificare, descompunere, deshidratare, maleabilizare, metal injection moulding (MIM), pyrolyses, rapid prototyping, recoacere, sintering, sinteză, sublimation

CARACTERISTICI STANDARD

- | 1600 and 1700 °C maximum operating temperatures
- | Sistem de control programabil, model EPC3016P1
- | Protecție la supra-încălzire
- | Designed for vertical use
- | Accepts work tubes with outer diameters up to 100 mm for use with modified atmosphere
- | Heated lengths of 120, 250, 500 or 700 mm
- | High grade type B thermocouple
- | Low thermal mass ceramic fibre insulation
- | Vertically hanging high quality MoSi₂ heating elements
- | Supplied with separate control box and 3 m cable, plug and socket
- | Conexiune Ethernet

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
DETALII TEHNICE

View inside the HTRV-A split tube furnace

1. Outer case
2. Ceramic fibre end insulation
3. Ceramic fibre case insulation
4. Heating elements
5. Ceramic fibre inner insulation

View inside the HTRV-A split tube furnace

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