



### Technical specification information for Carbolite CTF 1700 & 1800 tube furnaces

Model	CTF 17/-/300	CTF 17/-/600	CTF 18/-/300	CTF 18/-/600
Maximum temperature (°C)	1700	1700	1800	1800
Tube length (mm)	650	950	650	950
Heated length (mm)	300	600	300	600
Uniform length (mm - $\pm 5^{\circ}\text{C}$ )	200	400	200	400
Tube inner diameter (mm)	25/50/75	25/50/75	25/50/75	25/50/75
Power rating (kW)	6	9	6	10
Outer measurements excluding tube (mm - h x w x d)	755/600/555	900/950/630	755/600/555	755/900/555
Weight (kg)	126	220	130	230