



Technical specification information for Carbolite PIC cooled incubators

Model	PIC 30	PIC 60	PIC 120	PIC 200
Temperature range (°C)	0 - 60			
Maximum continuous working temperature (°C)	60			
Volume (l)	30	60	120	200
Chamber dimensions (mm - h x w x d)*	300 x 290 x 320	400 x 392 x 420	500 x 492 x 520	750 x 592 x 520
External dimensions (mm - h x w x d)	715 x 705 x 510	815 x 805 x 570	915 x 905 x 670	1165 x 1005 x 670
Number of shelves supplied	2	2	2	2
Temperature uniformity @ 37°C (°C)	± 1			
Temperature stability (°C)	± 0.25			
Power rating at 240V (W)	250	700	700	1000