



Technical specification information for Carbolite RHF 1400°C – 1600°C range of furnaces

Model	Maximum operating temperature (°C)	Volume (litres)	Internal dimensions (mm - h x w x d)	External dimensions (mm – h x w x d) with door closed	Heat up time to 100°C below maximum temperature (mins)	Maximum power (kW)	Holding power (kW)
RHF 14/3	1400	2.9	120 x 120 x 200	655 x 435 x 610	33	4.5	1.9
RHF 14/8	1400	7.8	170 x 170 x 270	705 x 505 x 675	22	8.0	3.2
RHF 14/15	1400	14.8	220 x 220 x 305	810 x 690 x 780	35	10.0	2.9
RHF 14/35	1400	35	250 x 300 x 465	885 x 780 x 945	38	16.0	6.0
RHF 15/3	1500	2.9	120 x 120 x 200	655 x 435 x 610	45	4.5	2.0
RHF 15/8	1500	7.8	170 x 170 x 270	705 x 505 x 675	40	8.0	3.5
RHF 15/15	1500	14.8	220 x 220 x 305	810 x 690 x 780	45	10.0	3.0
RHF 15/35	1500	35	250 x 300 x 465	885 x 780 x 945	46	16.0	6.2
RHF 16/3	1600	2.9	120 x 120 x 200	655 x 435 x 610	55	4.5	2.3
RHF 16/8	1600	7.8	170 x 170 x 270	705 x 505 x 675	60	8.0	4.0
RHF 16/15	1600	14.8	220 x 220 x 305	810 x 690 x 780	60	10.0	3.5
RHF 16/35	1600	35	250 x 300 x 465	1530 x 900 x 1020	90	16.0	5.0